

smartPillCap Instructions

Features

Comes with 3 pre-programmed options:

8 Hours

12 Hours

24 hours.

General operation

The smartPillCap turns on when the alarm is set for a time interval.

When the alarm goes off, it can be cancelled by pressing the button for less than 5 seconds. The alarm will always go off at the time it was set for, even if it was stopped later. For example, if a 24 hour alarm goes off at 9AM on Monday and is stopped at 9:05AM; with the short press of a button, this same alarm will sound again on Tuesday morning at 9:00AM.

The alarm on this device can be set to any start time. Simply press it at the time you want to be reminded.

To set the alarm for 8AM, press the button on top of the device and hold it for 5 seconds at 8 AM.

If you want to change it to 9AM, hold the button for 5 seconds at 9AM.

When the alarm goes off, it will follow the profile until the next alarm time.

Button operation

There are three button press types:

Short Press: Used to stop an alarm condition until the next alarm time.

Less than 5 seconds will be answered with 2 short beeps.

Medium Press: Used to synchronize the alarm

Hold the button for 5 seconds and you will hear a long beep. This is to let you know that the alarm has been synchronized.

Long Press: Used to change alarm time.

The button is held for 10 seconds. There will be a long beep at 5 seconds and two beeps at 10 seconds. Release the button to set the new alarm time. The device will beep twice when the time has been set.

Alarm Operation

Level 1 alarm -This is your signal. You will hear three beeps every 3-4 seconds. The LED light flashes on the same interval until turned off.

Level 2 - The device flashes every 2 minutes. It also beeps every 5 minutes for a total of 2 hours until turned off.

Level 3 - LED flashes every 5 minutes and the device beeps every 10 minutes, until the next alarm time or until turned off.